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de

Instruction manual

en

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testo 230
pH- and temperature measuring instrument

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de

Instruction manual

en


| | |
|---|----|
| Instrument description | 3 |
| Display/keyboard..... | 4 |
| Operation instructions | 5 |
| Safety instructions | 5 |
| Connection variations | 6 |
| pH/°C measurement..... | 6 |
| pH measurement..... | 7 |
| °C measurement..... | 7 |
| Redox measurement..... | 7 |
| Instructions of use | 8 |
| Connection of probes..... | 8 |
| Switching on the instrument..... | 9 |
| Auto-off function..... | 9 |
| 1 point calibration | 10 |
| 2 point calibration | 12 |
| pH measurement | 15 |
| Hold function | 16 |
| Holding measured values..... | 16 |
| Auto-hold function..... | 16 |
| Configuration menu | 17 |
| °C/°F conversion • Auto off function..... | 17 |
| Clear function..... | 18 |
| Care and maintenance | 19 |
| Battery replacement..... | 19 |
| General remarks..... | 20 |
| Error messages | 21 |
| Error 1-3..... | 21 |
| Error 4-6..... | 22 |
| Correcting measuring errors | 23 |
| Unstable display..... | 23 |
| Low slope and/or slow response time..... | 23 |
| Drifting potential..... | 23 |
| Technical data | 24 |
| Instrument..... | 24 |
| pH electrodes/temperature probe..... | 25 |
| Warranty..... | 25 |
| Ordering data | 26 |
| Appendix | 27 |
| Which pH electrode for which application..... | 27 |
| pH scale..... | 28 |
| Buffer table..... | 28 |

Instrument description

Display /keyboard

Display of measured value (pH / mV)
(°C for temperature measurements without
pH electrode)

Unit of temperature measurement

Probe recognition

- MANelectrodes without temp. sensor
- Symbol + AUTO...external temp. sensor
- AUTOelectrode with temp. sensor

DIN/NBS buffer (display in calibration menu)


Battery control

When switching on: battery voltage
During operation: low battery warning

HOLD function/Auto-hold function

Display of calibration menu

Auto-off display

Operation instructions

testo 230

measures pH values and redox voltages in aqueous solutions and in semi-solids to solids, as well as temperature.

IP 54

testo 230

is a simple to clean, rugged and easy-to-use pH and temperature measuring instrument covering a range of pH 0 to 14, 0 to 2000 mV as well as -50 to +150 ° C (-60 to 300 °F).

The water-resistant housing (IP 54 protection, showerproof and dust protection), its easy use and the range of different electrodes (probes) makes the instrument suitable for a wide range of applications.

Safety instructions

Please read before using instrument!

Do not measure on live parts!

Please read the "**testo pH electrodes**" instruction manual before using the instrument!

Observe permissible storage, transport and operating temperatures (e.g. protect measuring instrument from direct sunlight).

IP 54 protection is only guaranteed when probes or filler plugs are connected.

Electrolyte solutions are corrosive, do not bring into contact with the measuring instrument or other sensitive surfaces !

Ether, ester, ketone along with aromatic or halogenated hydrocarbons attack the shaft material and should, therefore, not come into contact with the electrode.


The electrode housing is made partially of glass,
ATTENTION, RISK OF BREAKAGE!

Should the instrument be opened, improperly handled or if force is applied, the warranty will no longer be valid.

Connection variations

pH/°C measurement


The pH value and temperature can be measured with the type 04 pH and type 05 pH electrodes. The temperature of the solution is automatically incorporated in the measurement and calibration procedure.

04 pH

05 pH

The pH value can be measured with the type 01 pH, 02 pH, 03 pH and 13 pH electrodes while the temperature of the solution can be measured with the type 02T, 24 T, 70T, 22T and 09T NTC probes. The temperature of the solution is automatically incorporated in the measurement and calibration procedure.

01 pH

02 pH

03 pH

13 pH

02 T

09 T

22 T

70 T


24 T

Connection variations

pH measurement

The pH value can be measured with the type 01 pH, 02 p, 03 pH and 13 pH electrodes. The temperature of the solution must be entered manually during measurement and calibration.


01 pH

02 pH

03 pH

13 pH

°C measurement

By connecting the type 02T, 24 T, 70T, 22T and 09T NTC probes, the testo 230 becomes a fully-competent thermometer covering the range -50 °C to +150 °C.

02 T

09 T


22 T

70T

24 T

Redox measurement

The redox voltage of the solution can be measured by connecting the type 06 mV redox electrode. An NTC probe can be used to measure temperature, although it has no effect on the redox measurement.

06 mV

Instructions of use


Connection of probes

The cable from the electrode must be connected before the instrument is switched on.

Insert the plug into the socket and turn to the right.


When connecting an external probe, pay attention to the markings (><) on the instrument and the probe plug.

Never pull on the lead when removing the probe!


The electrodes or temperature probe can be used loose or attached to the side of the instrument.

Should no external probe be used, filler plugs should be used to protect the socket and the measuring instrument.

Instructions of use

Switching on the instrument


Switch on instrument.

A short display and segment test then follows, after which a battery control test indicates the battery voltage (see also battery replacement).

The measuring instrument is now operational.

ATTENTION!
Calibration must be carried out before the first measurement (and then in intervals appropriate to the application)!

Auto-off function


Auto-off function

If the instrument is not operated for 10 minutes, the instrument switches off automatically

The auto-off function appears in the display. The circle disappears from the display at a rate of 2.5 minutes/segment. When the last segment disappears, the instrument switches itself off. (Switch off Auto-off function, see chapter on "Configuration menu").

1 point calibration


The 1 point calibration enables the instrument to adapt quickly to the electrode. Calibration (1 point or 2 point) must be carried out before the first measurement.

The following buffers are, thus, available:
buffer pH 2; pH 7, pH 10
DIN buffer pH 1.68; pH 4.008*;
pH 6.865*; pH 9.18


This procedure should be repeated before every measurement when measurements are to be carried out in extreme solutions, when the electrode is used frequently or when a very high accuracy is required.

Connect the electrode to the instrument. Open the vent plugs (if present) on the electrode. Carefully remove the immersion cap from the electrode. Rinse the electrode in water, dry it and place it in the buffer solution.

Switch on the instrument.


Select the "CAL" menu with the page key. "CAL" flashes in the display.

Confirm with **SELECT**.

The temperature is displayed.

Should no external probe be connected and should an electrode without temperature sensor be used, the temperature must be entered manually with the arrow keys (▲▼). Wait 5 seconds, the instrument then accepts this temperature value.

*Values displayed: pH 4.01; pH 6.87

1 point calibration

The instrument now asks for the buffer. Select the buffer required with the arrow key (▲) and start the calibration procedure with **SELECT**. The following solutions can be used: pH 1.68 DIN/NBS; 2.00; 4.00; 4.008* DIN/NBS; 6.865* DIN/NBS; 7.00; 9.18 DIN/NBS; 10.00

"CAL" flashes during calibration (the calibration lasts for min. 20 sec. and max. 10 min.). A signal sounds the end of the calibration procedure. Should the calibration procedure end abruptly after max. 10 minutes, the electrode is probably faulty. Please contact our service department.

The instrument asks for a second buffer. This menu is part of the 2 point calibration and is ignored here.

Leave the menu by pressing the page key. The set slope value is displayed "-58.0mV/pH".

Confirm with **SELECT**.


The asymmetry voltage is displayed. Confirm with **SELECT**, the instrument then switches to the pH measuring mode.

The 1 point calibration is now completed.

You can now begin the measurement.

*Values displayed: 4.01 DIN/NBS; 6.87 DIN/NBS

2 point calibration

The 2 point calibration enables the instrument to adapt accurately to the electrode. Calibration (1 point or 2 point) must be carried out before the first measurement.

The following sets are available for calibration:
buffer set pH 2.00/pH 4.00
buffer set pH 7.00/pH 10.00
DIN-buffer pH 1.68; pH 4.008*;
pH 6.865*; pH 9.18

The difference between buffer solutions A (1 bu) and B (2 bu) must be greater than 0.32 pH.

This procedure should be repeated before every measurement when measurements are carried out in extreme solutions, when the electrode is used frequently or when a very high accuracy is required.

Connect the electrode to the instrument.

With electrodes Type 02 and 05:

Open the vent plugs on the electrode. Carefully remove the immersion cap from the electrode. Rinse the electrode in water, dry it and place it in the first buffer solution (A).



Switch on the instrument.

Select the "CAL" menu with the page key. "CAL" flashes in the display.

Confirm with **SELECT**.


The temperature is displayed.

Should no external probe be connected and should an electrode without temperature sensor be used, the temperature must be entered manually with the arrow keys (▲▼). Wait 5 seconds, the instrument then accepts this temperature value.

*Values displayed: pH 4.01; pH 6.87

2 point calibration


The instrument now asks for the first buffer. Select the buffer required with the arrow key (▲) and start the calibration procedure with **SELECT**. The following solutions can be used:
pH 1.68 DIN/NBS; 2.00; 4.00; 4.008* DIN/NBS; 6.865* DIN/NBS; 7.00; 9.18DIN/NBS; 10.00

"CAL" flashes during calibration (the calibration lasts for min. 20 sec. and max.10 min.). A signal sounds the end of the calibration procedure. Should the calibration procedure end abruptly after max. 10 minutes, the electrode is probably faulty. Please contact our service department.

Once this procedure is completed the electrode is rinsed in water, dried and then placed in the second buffer solution (B).


Should the required buffer not correspond with the set value, reset the value with the arrow key (▲).

Start the calibration procedure with **SELECT**.

*Values displayed: 4.01 DIN/NBS; 6.87 DIN/NBS

2 point calibration

"CAL" flashes during calibration (the calibration lasts for min. 20 sec. and max. 10 min.). A signal sounds the end of the calibration procedure. Should the calibration procedure end abruptly after max. 10 minutes, the electrode is probably faulty. Please contact our service department.


The slope value is displayed.


Confirm with **SELECT**.

The asymmetry voltage is displayed. Confirm with **SELECT**, the instrument then switches to the pH measuring mode.

The 2 point calibration is now completed.

You can continue with the measurement.

Connect the electrode to the measuring instrument.

Open the vent plugs on the electrode (for pressure compensation).

Remove the immersion cap from the electrode.

Rinse the electrode in water and dry off (see also pH electrodes instruction manual).


Place the electrode (and the external temp. probe) in the measuring solution.

Switch on instrument.


The display indicates the pH value (at first unstable) and temperature.

Should no external probe be connected and should an electrode without temperature sensor be used, the temperature must be corrected manually with the **arrow keys**. First activate the **SELECT key**, "MAN" flashes in the display, then correct the value with the **arrow keys** (▲▼). Wait 5 seconds, the instrument then accepts this temperature value.

The value displayed is correct when it has stabilised.

Press the **page** key to enter the following menu points:

- **Auto-Hold**, see chapter on hold function.
- **Redox measurement**, for measuring redox voltage in conjunction with the type 06 mV redox electrode.
The electrode voltage is displayed when a pH electrode is connected.
- **Temperature measurement**, the instrument becomes a fully-competent temperature measuring instrument when a NTC probe is connected.
- **Calibration menu**, see 1 / 2 point calibration.

Press the **page** key to return to the measuring menu.

Hold function

Holding measured values

The **HOLD** key can be used to hold a measured value in display.

Press the **HOLD** key.

"Hold" appears next to the measured value. The measured value can be held in display as long as necessary, the Auto-off function is deactivated.

Press the **HOLD** key again to end this function, the current measured value is displayed.

Auto-Hold function

The Auto-Hold function monitors the response time during a measurement. An optical and acoustic signal indicates that the measured value has stabilised.

Start the measuring procedure (see pH measurement).

Select the "Auto-Hold" with the page key.

"Hold" flashes in the display. A signal sounds to indicate that the measured value has stabilised, the measured value and "Hold" are held in display. The measuring procedure is completed.


Press the **HOLD** key to return to the normal measuring mode.

Should the response time be > 10 min, "ERR 6" is displayed (see error messages).

Configuration menu

°C / °F conversion • Switching auto-off function on or off

testo 230 enables the conversion of the unit of measurement from °C to °F and vice-versa. The auto-off function can also be switched off.

The auto-off function increases the battery life and thus protects the environment. Should the auto-off function be disconnected, do not forget to switch the instrument off!

Keep the **HOLD** key pressed when switching on the instrument and press the **I/O** key once.

Keep the HOLD key pressed until the segment test has been completed.

"°C" or "°F" flashes in the display.

Set the required unit of measurement with the **SELECT** key and confirm with **HOLD**.

The "Auto-off" symbol flashes in the display.


SELECT can be used to switch this function on or off.

Confirm the function with **HOLD**.

The "Clear" function now appears in the display.

Configuration menu

Clear function

The "Clear" function can be reached after the °C/°F and the Auto-Off modes.

In this menu, the values stored during calibration can be altered, for example, when you have to replace an electrode and calibration is not possible.

In this menu, the slope and asymmetric voltage values are set to the ideal values of -58 mV/pH and 0 mV and stored.

Call up "clr".

Activate "yes" with the **SELECT** key.

Leave the menu with **HOLD**.

The previously-stored values are erased and the above-mentioned ideal values are stored.

The instrument switches to the measuring menu.

Should "no" be selected in the clear mode and confirmed, the previously-stored values are maintained.

Care and maintenance

Battery replacement


The battery voltage is displayed after the segment test each time the instrument is switched on and can, thus, be checked regularly.

Should the battery voltage drop below the permissible value, "BAT" flashes in the display.

Should the battery voltage no longer be sufficient to guarantee the correct functioning of the instrument, the measuring instrument switches off automatically.

Open the battery compartment.

Remove the old battery and insert the new battery, type 9V (IEC 6 F 22)

When replacing the battery, take care not to damage the ring.

(observe polarization).

Care and maintenance

General remarks

testo 230 requires no maintenance and is therefore not subject to particular maintenance intervals.

Clean the housing with a damp cloth. Weak household cleaners can be used.

The instrument is suitable for contact with weak acids and electrolyte.

Strong detergents or solvents should never be used for cleaning the measuring instrument.

After use in an aggressive environment, the instrument should be cleaned with a damp cloth.

Please note that the IP 54 protection is only guaranteed when original plugs or filler plugs are connected.
Never clean with a water jet.

IP 54

Read the instruction manual "**testo pH electrodes**" for the care and maintenance of pH electrodes.

pH or mV meas. range has been exceeded

The pH measured value is not within the range 0 to 14 pH.

- Possible causes and • remedies:
 - no pH electrode is connected.
 - connect pH electrode.
 - incorrect calibration.
 - recalibration.
 - faulty pH electrode.
 - check pH electrode and replace if necessary.
 - lower pH value with electrode with high slope and asymmetry voltage and simultaneously high measuring temperature.
 - reduce measuring temperature.


The mV measured value is not within the range 1999 to +1999 mV.

- Possible causes and • remedies:
 - no electrode is connected.
 - connect electrode.
 - faulty redox electrode
 - check electrode and replace if necessary.

Temp. range of the instrument has been exceeded.

- Possible causes and • remedies:
 - temperature < -50 °C or > +150 °C.
 - remove probe from substance being measured.

Temp. range has been exceeded during calibration.

- Possible causes and • remedies:
 - temperature of the buffer solution is outside the permissible range.
 - interrupt calibration by pressing the PAGE key (the values from the last calibration will remain stored).
Repeat calibration with suitable temperature.

Slope of the electrode is not within the range -50 mV/pH to -62 mV/pH

- Possible causes and • remedies:
 - faulty electrode.
 - interrupt calibration by pressing the PAGE key (the values from the last calibration will remain stored).
- Use new electrode.


Zero point of the electrode is not within the range -60 to +60 mV.

- Possible causes and • remedies:
 - faulty electrode.
 - interrupt calibration by pressing the PAGE key (the values from the last calibration will remain stored).
- Use new electrode.

Response time of the electrode too long
(during Auto-Hold or calibration).

- Possible causes and • remedies:
 - electrode is dirty
 - electrode faulty.
 - erase error message by pressing the PAGE key (interrupt calibration by pressing the PAGE key, the values from the last calibration will remain stored).
- rinse the electrode and continue the ,
procedure, should the error reoccur, use a
new electrode.

Internal instrument error.

- Possible causes and • remedies:
 - faulty measuring instrument.
 - return instrument to Testo's service department.

Correcting measuring errors

Unstable display

| Cause | Remedy |
|--|---|
| - Air bubbles in glass ball | • Knock carefully on the electrode with your fingers or swing the electrode round in a circle. |
| - Air bubbles in KCl crystals | • See above. • Heat electrode to 60 °C in water bath until sufficient KCl is dissolved and the air bubbles have disappeared. |
| - Bad contact | • Check the cable connection. |
| - Low ionic strength of test solution (e.g. part-aqueous solution) | • Use the type 05 pH electrode. |

Low slope and/or slow response time

| Cause | Remedy |
|--|--|
| - Calibration error caused by old buffer | • Use fresh buffer solutions of a high quality for calibration. |
| - Contaminated glass membrane | • Clean the glass ball with soft paper or lense cloth and leave the electrode membrane to soak in a pH 4.0 buffer. • Clean the electrode in tap water and rinse with distilled water. Then leave the electrode in a pH 4.0 buffer solution with an open KCl vent plug for more than 30 minutes, stirring from time to time. |

Drifting potential

| Cause | Remedy |
|--|--|
| - Oil and grease contaminated. | • Clean the electrode with a water-mixable solution (e.g. acetone) and rinse the electrode in warm, distilled water. |
| - Sulphide contamination of porous pin | • Leave the electrode to soak for 24 hours in a thiourea solution of 0.1 M HCl and rinse with distilled water. |

Technical data

Instrument

| | |
|---------------------------------|--|
| Measuring range: | pH 0 to 14 -50 to +150 °C ±1999 mV |
| Resolution: | pH 0.01 0.1 °C (°F) 1 mV |
| Accuracy: ± 1 digit | ±0.01 pH ±0.4 °C/0.7°F (-50 to -25 °C/-60 to -10 °F) ±0.2 °C/0.4°F (-25 to +75 °C/-10 to 165 °F) ±0.4 °C/0.7°F (+75 to +100 °C/+165 to +210 °F) ±0.5 % of m.v. (+100 to +150 °C/+210 to +300 °F) ±1 mV (0 to ±999 mV) ±2 mV (±999 to ±1999 mV) |
| Compensation de temp.: | man -10 to +150 °C (-15 to +300 °F) auto -50 to +150 °C (-60 to +300 °F) |
| Operating temp.: | 0 to +40°C (+32 to +105 °F) |
| Storage/transport temp.: | -20 to +70°C(-4 to +160 °F) |
| Display: | two lines |
| Connections: | <ul style="list-style-type: none">• BNC compatible socket for pH or pH/°C• Mini DIN socket for temp. probes (testo 110) |
| Battery life: | approx. 100 hours. |
| Weight: | 180 g (incl. battery) |
| Other: | IP 54 low battery warning auto switch off (°C - °F conversion) |
| Housing material: | ABS |

Technical data

Electrodes/temperature probes

| Type | no. | Application description | | |
|------------------------------------|-----|---|---|-----------------------------------|
| pH with- out °C | 01 | Universal electrode for general liquid solutions | - pH 0 to 14 - plastic shaft - fixed cable with BNC | 0 to 60 °C (short-term 80°C)* |
| | 02 | Laboratory electrode , a chemical-resistant and highly-accurate electrode | - pH 0 to 14 - glass shaft - plug-in head**) | 0 to 80 °C (short-term 100°C)* |
| | 03 | Penetration electrode for semi-solids such as yoghurt, jam, milk | - pH 2 to 14 - glass shaft - spec. polymer | 0 to 40 °C (short-term 60°C)* |
| | 13 | Robust penetration probe for solids and semi-solids such as meat, fruit, bread | - pH 2 to 14 - glass tip with plastic shaft - spec. polymer | 0 to 40 °C (short-term 60°C)* |
| pH with °C | 04 | Universal electrode with integrated temperature sensor - accurate pH values are obtained automatically! | - pH 0 to 14 - spec. plastic shaft - fixed BNC compatible cable - temperature sensor | 0 to 60 °C |
| | 05 | Special electrode probe for atmospheres low in ions (drinking water) with integrated temp. sensor | - pH 0 to 14 - single pore diaphragm - fixed BNC compatible cable | 0 to 60 °C |
| Redox | 06 | Redox electrode The special electrolyte link supplies reproducible meas. values in environmental meas. technology | - ± 1999 mV - glass shaft - plug-in head**) | 0 to 80 °C |
| Temperature | 02 | Immersion/penetration probe | - -50...+150 °C | - t ₉₉ = 6 sec. |
| | 99 | Surface probe | - -50...+150 °C | - t ₉₉ = 25 sec. |
| | 22 | Food probe | - -50...+150 °C | - t ₉₉ = 10 sec. |
| | 70 | Laboratory probe with glass tube | - -50...+150 °C | - t ₉₉ = 40 sec. |
| | 24 | Robust penetration probe | - -50...+150 °C | - t ₉₉ = 10 sec. |

*) Short-term < 10 minutes

**) Plug-in electrodes require an electrode cable, part no. 0554.2317

Warranty

testo 230 measuring instrument.....12 months
Electrodes.....12 months*
Probes.....12 months

If the instrument is opened, improperly handled or if force is applied, the warranty will no longer be valid!

* for manufacturing faults

| Instrument | Part no. |
|--|------------------|
| testo 230 | 0560.2304 |
| testo 230 , Set 1 universal | 0563.2305 |
| testo 230 , Set 2 food..... | 0563.2306 |
| testo 230 , Set 3 low-cost | 0563.2307 |

Electrodes and probe

| | | | |
|---------------------------|---------|---------------------|------------------|
| Universal electrode | type 01 | pH without °C | 0650.0623 |
| Laboratory electrode | type 02 | pH without °C | 0650.1223 |
| Penetration electrode | type 03 | pH without °C | 0650.0225 |
| Penetration electrode | type 13 | pH without °C | 0650.0245 |
| Universal electrode | type 04 | pH with °C | 0650.0823 |
| Special electrode | type 05 | pH with °C | 0650.1623 |
| Redox electrode | type 06 | Redox..... | 0650.2523 |
| Immersion/penetration pr. | type 02 | °C | 0613.0211 |
| Surface probe | type 99 | °C | 0613.9911 |
| Food probe | type 22 | °C | 0613.2211 |
| Laboratory probe | type 70 | °C | 0613.7011 |

Accessories

| | |
|---|------------------|
| Electrode cable for plug-in electrodes S7-BNC | 0554.2317 |
|---|------------------|

Testo buffer sets

| | |
|---------------------------------------|------------------|
| pH 4/7/10 (50 ml each)..... | 0554.2320 |
| pH 4/7 (50 ml each)..... | 0554.2321 |
| pH 2 (50 ml)..... | 0554.2322 |
| pH 4.00; pH 7.00 (500 ml each) | 0554.2337 |
| pH 7.00; pH 10.00 (500 ml each) | 0554.2338 |

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| | |
|-----------------------------------|------------------|
| pH 4.008/6.865 (250 ml each)..... | 0554.2339 |
| pH 6.865/9.180 (250 ml each)..... | 0554.2340 |

Storage and refill solution

| | |
|---|------------------|
| Storage solution; (50 ml) | |
| for type 01, 02, 04, 06 (refill solution for type 02) | 0554.2332 |
| for type 03/13 electrode..... | 0554.2318 |
| Storage and refill solution; (50 ml) | |
| for type 05 electrode..... | 0554.2319 |
| Redox standard solution 358 mV; 50 ml | 0554.2333 |

Which electrode for which application


Key: - not suitable, o conditionally suitable, + suitable

| Type | 01 | 02 | 03 | 04 | 05 | 13 |
|---|----|----|----|----|----|----|
| Application | | | | | | |
| Sewage samples | + | o | | + | + | |
| General aqueous substances | + | + | | + | + | |
| Aquariums | + | + | | + | o | |
| Beer, fruit juices, wine | o | o | | o | + | |
| Yoghurt, cheese | | | + | | | + |
| Substances containing egg white | - | - | + | - | o | + |
| Aqueous emulsions | | | | | + | |
| Part-aqueous emulsions | | | | | + | |
| Earth (suspensions) | o | o | | o | + | |
| Extreme pH values (pH<1, >13) | o | o | - | o | + | |
| Penetration meas. in meat | - | - | o | - | - | + |
| Penetration meas. in fruit, vegetables | - | - | o | - | - | + |
| Substances containing hydrofluoric acid | - | o | - | - | o | |
| Galvanic sewage | | | | | + | |
| Hot electrolyte | - | o | - | - | - | |
| Highly-viscous solutions | | o | o | | - | |
| Infusion solutions | | o | | | + | |
| Solutions low in ions | - | o | | o | o | |
| Jams | | | + | | - | + |
| Cosmetic products | | o | | | - | |
| Leather, paper | o | o | | o | + | |
| Milk | | o | + | | | + |
| Rain water | - | o | | o | o | |
| Cream | | o | + | | | + |
| Brine | o | + | | + | + | |
| Swimming pools | + | + | | + | + | |
| Soaps, detergents | | | | | + | |
| Aqueous suspensions | | | | | + | |
| Part-aqueous suspensions | | | | | + | |
| Pastry, bread | | | + | | - | + |
| Part-aqueous suspensions, > 10 % H ₂ O | | o | | | + | |
| Part-aqueous suspensions, < 10 % H ₂ O | | | | | o | |
| pH measurement for temp. up to 100 °C | - | + | - | - | - | |
| Tris buffer solutions | o | o | - | o | + | |

Special electrodes are available on request for special applications: lacquer and paint, highly-viscous solutions, photographic solutions, measurements on pure surfaces and for measurements at T < 0 °C.

Appendix

pH scale

Buffer table

The temperature-dependant changes in the buffer solution are listed in the table below.
 (DIN buffer according to DIN 19266)
 This is automatically taken into consideration during calibration.

| °C | DIN Buffer | | | | Testo Buffer | | | |
|----|------------|-------|-------|-------|--------------|------|------|-------|
| | 1.685 | 4.006 | 6.865 | 9,180 | 2.00 | 4.00 | 7.00 | 10.00 |
| 10 | 1.670 | 4.000 | 6.923 | 9,332 | 2.00 | 4.00 | 7.07 | 10.18 |
| 15 | 1.672 | 3.999 | 6.900 | 9,276 | 2.00 | 4.00 | 7.04 | 10.14 |
| 20 | 1.675 | 4.001 | 6.881 | 9,225 | 2.00 | 4.00 | 7.02 | 10.06 |
| 25 | 1.679 | 4.006 | 6.865 | 9,180 | 2.00 | 4.00 | 7.00 | 10.00 |
| 30 | 1.683 | 4.012 | 6.853 | 9,139 | 1.98 | 4.01 | 6.99 | 9.95 |
| 35 | 1.688 | 4.021 | 6.844 | 9,102 | 1.99 | 4.02 | 6.98 | 9.91 |
| 38 | 1.691 | 4.027 | 6.840 | 9,081 | | | | |
| 40 | 1.694 | 4.031 | 6.838 | 9,068 | 2.00 | 4.03 | 6.97 | 9.85 |
| 45 | 1.700 | 4.043 | 6.834 | 9,038 | | | | |
| 50 | 1.707 | 4.057 | 6.833 | 9,011 | 1.99 | 4.05 | 6.96 | 9.78 |
| 55 | 1.715 | 4.071 | 6.834 | 8,985 | | | | |
| 60 | 1.723 | 4.087 | 6.836 | 8,962 | 1.99 | 4.08 | 6.96 | 9.75 |

testo AG

Postfach 11 40, 79849 Lenzkirch

Testo-Straße 1, 79853 Lenzkirch

Telefon: (07653) 681-0

Fax: (07653) 681-100

E-Mail: info@testo.de

Internet: <http://www.testo.com>