Installation
Installation
Installation

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Table of Contents

Introduction	5
Survey	5
System requirements	5
Hardware requirements	5
Software requirements	5
Unpacking/Scope of delivery	7
Installation of Software	9
Putting in CD-ROM	11
Starting Setup assistant	11
Give user group	13
Choose installation index	12
Start installation	15
Connection of Interface (write-/reading device)	17
General	17
Connecting interface to computer	19
Solving installation problems	21
First steps with „Winlog.med Validation“	25
Survey	25
Adjusting programme	27
Preparation of validation	27
Evaluation of validation	29
Looking into validation report	31

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Introduction

Survey

Winlog.med Validation is a reporting and evaluation software. It meets the strict requirements of pharmaceutical validation, qualification and measuring tasks. The software has a programme surface which is easy to operate and offers numerous possibilities of generating meaningful reports.

System requirements

So that the software can be operated on your computer without any problems, your computer has to fulfil the following requirements.

Hardware requirements

- Processor speed at least 1,5 GHz
- 512 MB working storage
- 200 MB free hard disk storage
- USB (Universal Serial Bus)

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Unpacking/Scope of delivery

- ◆ Check the packaging of the equipment and the contents after reception to make sure they are undamaged. Also check whether the contents of the equipment packaging corresponds to your order.

Scope of delivery:

These installation-instructions

CD-ROM

Interface

Interface-Netadapter

USB-cable

If you have any reason for a complaint, please contact us. Find our contact data on the reverse of the envelope.

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Installation of Software

Winlog.med Validation has its own installation programme, the setup assistant, which guides you in several steps through the installation of the software. Therefore the installation of the software is very easy:

- put in CD-ROM
- start Setup-Assistant
- give user group
- choose installation index
- start installation

In the following the separate installation steps are briefly described.

In the case of problems the setup assistant gives you corresponding hints with suggestion solutions. Therefore in every case observe the messages of the setup assistant on the monitor!

After the installation of the software you just have to connect the interface (write-reading device) (see page 16).

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Putting in CD-ROM

- ◆ Put the CD-ROM into the respective drive of your computer.

Starting Setup-Assistent

If the autostart function is switched on on your computer, the setup assistant will start automatically after a few moments.

If the autostart function is switched off on your computer, you have to start the setup assistant manually.

- ◆ In this case select „Start: Execute „, from the Windows tool bar
- ◆ Now indicate the drive where the CD-ROM is (e.g. „D:“) and directly write after it „\setup.exe“.

In case it is necessary, you can cancel the installation any time. For this click on the button „Cancel“ and confirm the query of the setup assistant with „Yes“.

- ◆ In order to continue with the installation, click on the button „Next“.

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
After clicking on „OK“ the Setup-Assistent starts. Now the setup assistant guides you through the installation of the software.

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Give user group

You can install the software in such a way that only you have an access to the software. If you are system administrator, you can also indicate in installation, that all users of the computer are permitted to have access to the software.

After this decision, the setup assistant asks you

Presetting: All users

Choose installation index

You can install the software in any index you want. The plant has already asser-
ted the index

C:\Programme\Winlog.med Validation

◆ **In order to choose another, already existing index**, click on „Browse“ and then choose the desired index.

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- ◆ **In order to indicate a not yet existing index and have it created by the setup assistant**, enter the index to be created manually into the line „Directory“.

Now all necessary preparations have been met.

- ◆ In order to start the installation, click on „Next“.

Start installation

When all installation requirements have been fulfilled, the setup assistant now installs the software completely.

In the case of problems, the setup assistant gives you corresponding hints with suggestion solutions. Therefore in every case observe the messages of the setup assistant on the monitor!

If everything has functioned without problems, then after some time the setup assistant will give the message, that the software was successfully installed and the installation is now completed.

After the installation of the software you only have to connect the interface (write-reading device) (see next page).

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Connection of Interface (Write-/reading device)

General

The software „Winlog.med Validation“ creates validation plans and evaluates validations. Data loggers are aids for this, which are programmed by the software and after validation are selected from the software.

Write-/reading devices (Interfaces) are needed for programming and to select the data loggers. These are connected over the USB (Universal Serial Bus) with the computer. This can also happen during the operation of the computer without the computer operation being disturbed because of this.

The operating system automatically recognises whether there are modifications to the Universal Serial Bus (USB) and gives respective messages on the monitor.

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
Connecting Interface to computer

Write-/reading devices (Interfaces) are connected with the computer over the USB (Universal Serial Bus). In order to connect your interface with the computer, proceed as follows:

- ◆ Ensure that the installation CD is in the computer drive. This contains control programmes (so called „drivers“) which are necessary for the operation of the interfaces.
- ◆ In case you are using an eightfold interface, connect this over the included net adapter to the power supply and switch it on (switch lights up).
- ◆ Connect the interface with the help of the included USB cable to the computer.

The operating system automatically recognises whether there are modifications to the USB and gives corresponding messages on the monitor screen.

The operating system treats the eightfold interface like eight separate hardware devices. The message „New hardware found“ therefore appears 16 times subsequently for the eightfold interface.

Shortly after the message „New hardware found“ the „Assistant for searching for new hardware“ appears on the monitor screen automatically.

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
Now the „Assistant for searching for new hardware“ installs the driver.

◆ If the installation has been successful, click on „Complete“.

In case you have connected the eight-fold interface, the procedure „New hardware found“ runs through another 15 times up to „Complete“.

To solve installation problems

In case it should come to problems in the installation of the drivers, the adjacent message appears.

In this case advance in the following manner:

◆ Select **Start → Preferences → Control panel → Administration → Computer administration**

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The window „Computer administration“ will open.

- ◆ Click there on „Hardware -Manager“.

Now on the right side appears a list of the hardware devices contained in the computer. All hardware devices which pose problems for the operating system are marked with a yellow symbol. In the adjacent picture this is the case with one of the COM-Ports.

- ◆ If you also discover a COM port with a yellow symbol on your computer, double click on the respective COM port name.

The window „Characteristics of USB Serial Port“ will open.

- ◆ Click on the button „Install Driver anew“.

Now the „Hardware-Update-Assistant“ opens and offers to install the driver.

- ◆ Proceed as described with the „Assistant for searching for new hardware“ on pages 18/19. The now necessary steps are identical with those described there.

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

First steps with „Winlog.med Validation“

Survey

After the double click on the programme symbol, the programme main window will open. This has two main areas: the area for the function selection (Navigator) and the working space.

The procedures of the programmes are subdivided into three areas: „Validations“, „Preferences“ and „Reports“.

In the area „Validations“ validation plans can be drawn up, process data can be evaluated according to a validation procedure or validation reports can be opened to be looked at.

The area „Preferences“ offers entry possibilities for

- general information (e. g. index for the source database, exsistant interfaces [write-/reading devices] for the data loggers or physical units to be used),
- the security (e. g. password, notification email addresses or electronical signature),
- the user administration (who may use the programme and with which rights?) as well as
- the hardware devices used (e. g. type of devices, manufacturer, hardware devices programme).

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In the area „Reports“ the recorded programme tasks can be browsed, filtered and shown as a report according to various features under „Audit Trail“.

Adjusting Programme

After the installation the programme needs certain basic information, in order to work in accordance with the regulations. But also afterwards in the current use of the programme, it can be necessary to change certain preferences.

Adjusted can be

- general information as for example the index of the database, the interface (write/reading device), the units of measurements to be used or information on the users firm.
- Information on security such as passwords, protection from misuse, email addresses for automatic notifications or electronical signature;
- for users with access to the programme and their user rights;
- the hardware devices used (hardware device manufacturer, type of device etc.).

Please obtain further details from the included Help on the monitor.

Preparation of Validation

The programme helps you to prepare a validation resp. plan one. It guides you to this in six subsequent masks (steps)

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through the preparation. In these masks it demands all necessary details to be given one after the other.

After having finished the details you can have the validation plan shown and printed out.

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Evaluation of validation

The programme helps you with the evaluation of a validation. After expiry of a validation cycle, it confirms whether a validation was successful or not according to the validation plan and the data from the data loggers. For this the programme selects the data from the data loggers and judges this according to the details given from the validation plan.

Each validation only lets itself be evaluated once. A validation report is made from the validation plan by the evaluation. This is available any time after the evaluation and can be looked into.

The programme guides you in three subsequent masks (steps) through the evaluation.

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Looking into the validation report

With the help of a report over a validation, each validation also allows itself to be judged at a later time and for example in the case of a dispute to be proven (securing of quality).

The programme directly files all data and details which are created in preparing a validation or in the evaluation of a validation, into the corresponding database. In calling up a validation plan or a validation report, the programme picks out the necessary data from the database and puts them into a prefabricated report form..

Reports of validations already evaluated are saved in databases belonging to the programme secure from manipulation. They can be looked into any time. In doing so the data are not changed.

You can directly print the report or e.g. export it into Portable Document Format (PDF).

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