Instruction Manual

HI 93748 Manganese Low Range ISM



This in

Dear Customer,

Thank you for choosing a Hanna product. This manual will provide you with the necessary information for the correct operation of the meter. Please read it carefully before using the meter. If you need additional technical information, do not hesitate to e-mail us at tech@hannainst.com.

This instrument is in compliance with **C€** directives EN 50081-1 and EN 50082-1.

PRELIMINARY EXAMINATION

Remove the instrument from the packing material and examine it carefully to make sure that no damage has occurred during shipment. If there is any damage, notify your Dealer.

Each Ion Specific Meter is supplied complete with

- · 9V Battery
- · Two Sample Cuvets and Caps
- · One Transport Cap

Note: Conserve all packing material until the instrument has been observed to function correctly. Any defective item must be returned in its original packing.

WARRANTY

HI 93748 is warranted for two years against defects in workmanship and materials when used for its intended purpose and maintained according to instructions.

This warranty is limited to repair or replacement free of charge.

Damages due to accident, misuse, tampering or lack of prescribed maintenance are not covered.

If service is required, contact the dealer from whom you purchased the instrument. If under warranty, report the model number, date of purchase, serial number and the nature of the failure. If the repair is not covered by the warranty, you will be notified of the charges incurred. If the instrument is to be returned to Hanna Instruments, first obtain a Returned Goods Authorization Number from the Customer Service department and then send it with shipment costs prepaid. When shipping any instrument, make sure it is properly packaged for complete protection.

To validate your warranty, fill out and return the enclosed warranty card within 14 days from the date of purchase.

All rights are reserved. Reproduction in whole or in part is prohibited without the written consent of the copyright owner, Hanna Instruments Inc., Woonsocket, Rhode Island, 02895, USA.

Hanna Instruments reserves the right to modify the design, construction and appearance of its products without advance notice.

The HI 93748 meter measures the manganese content in water and wastewater in the 0 to 300 $\mu g/L$ (ppb) range.


The meter uses an exclusive positive-locking system to ensure that the cuvet is in the same place every time it is placed into the measurement cell.

The reagents are in liquid and powder form and are supplied in bottles and in packets. The amount of reagent is precisely dosed to ensure maximum repeatability.

Display codes aid the user in routine operations.

The meters have an auto-shut off feature that will turn the instrument off after 10 minutes of non-use.

SPECIFICATIONS



SPECIFICATIONS

Range 0 to 300 μ g/L Resolution 1 μ a/L

Accuracy $\pm 2 \mu g/L \pm 3\%$ of reading

Typical EMC $\pm l \mu g/L$ Deviation

Light Source Light Emitting Diode @ 555 nm

Method Adaptation of the PAN Method. The reaction between manganese and the reagents

causes an orange tint in the sample

Light Detector Silicon Photocell

Environment 0 to $50 \square C$ (32 to $122 \square F$);

max 95% RH non-condensing

Battery Type/Life 1 x 9 volt/40 hours Auto-Shut off After 10' of non-use

Dimensions 180 x 83 x 46 mm (7.1 x 3.3 x 1.8")

Weight 290 g (10 oz.).

REQUIRED REAGENTS

<u>Code</u>	<u>Description</u>	<u>Quantity</u>
HI 93748 A -0	Ascorbic acid	2 packets
HI 93748 B -0	Alkaline-cyanide sol.	0.40 mL
HI 93748 C -0	0.1% PAN indicator	2 mL
HI 93703-51	Dispersing Agent	4-6 drops (only
	when i	necessary, see note)

REAGENT SETS

HI 93748-01 Reagents for manganese LR (100 pcs of HI93748A-0, 20 mL of HI93748B-0, 100 mL of HI93748C-0, 20 mL of HI93703-51)

HI 93748-03 Reagents for manganese LR (300 pcs of HI93748A-0, 60 mL of HI93748B-0, 300 mL of HI93748C-0. 60 mL of HI93703-51).

DISPLAY CODE GUIDE

This indicates that the meter is in a ready state and zeroing can be performed.

Sampling in Progress. This prompt appears each time the meter is performing a measurement.

This indicates that the meter is in a zeroed state and measurement can be performed.

A zero reading was not taken. Insert a sample before adding reagent and press 7FRO

Under range. A blinking "0.00" indicates that the sample absorbs less light than the zero reference. Check the procedure and make sure you use the same cuvet for reference (zero) and measurement.

Over range. A flashing value higher than the maximum concentration readable (see specifications) indicates that the sample absorbs too much light, meaning that the concentration is too high. Dilute the sample.

Light over range. The cuvet is not inserted correctly and an excess ambient light is reaching the detector. If the cover is properly installed, then contact your dealer or the nearest Hanna Customer Service Center.

Light under range. The zero sample is too dark for proper zeroing. If this is not the case, contact your dealer or the nearest Hanna Customer Service Center.

The "V" indicates that the battery voltage is getting low and the battery needs to be replaced.

This indicates that the battery is dead and must be replaced.

Note: once this indication is displayed, the meter will lockup. Change the battery to restart.

OPERATIONAL GUIDE

MEASUREMENT PROCEDURE

- Turn the meter on by pressing ON/OFF.



Note: for best results perform your tests between $20-24 \square C$.

- · When the LCD displays "- · ", it is ready.
- Fill one cuvet up to 1.5 cm (¾") below the rim with 10 mL of deionized water.




• Fill a second cuvet up to 1.5 cm (34") below the rim with 10 mL of sample.




 Add the content of one packet of HI 93748A ascorbic acid to each cuvet, replace the caps and shake gently.



 Add 0.2 mL of the HI 93748B alkaline-cyanide reagent solution to each cuvet, replace the caps and shake gently.







 Press READ TIMED and the display will show the countdown prior to zeroing the blank. Alternatively wait for 2 minutes and press ZERO. In both cases "SIP" will appear during zeroing.



Wait for a few seconds and the display will show "-0.0-". Now the meter is zeroed and ready for measurement



Insert the second cuvet with the reacted sample into the instrument.



Press READ DIRECT. "SIP" will appear during measurement.



The instrument directly displays concentration in μg/L of manganese on the display.

Note: a temperature above 30 □ C may cause turbidity. In this case, before zeroing and taking readings, add 2-3 drops of Dispersing Agent (HI93703-51) to each cuvet and swirl until complete dissolution of the turbidity.

INTERFERENCES

Interference may be caused by:
Aluminum above 20 mg/L
Cadmium above 10 mg/L
Calcium above 200 mg/L as CaCO₃
Cobalt above 20 mg/L
Copper above 50 mg/L
Iron above 10 mg/L
Lead above 0.5 mg/L
Magnesium above 100 mg/L as CaCO₃
Nickel above 40 mg/L
Zinc above 15 mg/L

TIPS FOR AN ACCURATE MEASUREMENT

The instruction listed below should be carefully followed during testing to ensure best accuracy.

- · Do not touch the cuvet walls with hands.
- In order to maintain the same conditions during the zeroing and the measuring phases, it is necessary to close the cuvet to prevent any contamination.
- Do not let the test sample stand too long after reagent is added or accuracy will be lost.
- Whenever the cuvet is placed into the measurement cell, it must be completely free of fingerprints, oil or dirt.
 Wipe it thoroughly with HI 731318 or a lint-free cloth prior to insertion.
- It is important that the sample does not contain any debris. This would corrupt the readings.
- It is possible to take multiple readings in a row, but it is recommended that a zero reading be taken for each sample and that the same cuvet is used for zeroing and measurement.
- It is important to discard the sample immediately after the reading is taken because the glass might become permanently stained.
- Shaking the cuvet can generate bubbles in the sample, causing higher readings. To obtain accurate measurements, remove such bubbles by swirling or by gently tapping the vial.
- All the reaction times reported in this manual are referred to 20□ C (68□ F). As a general rule of thumb, they should be doubled at 10□ C (50□ F) and halved at 30□ C (86□ F).

BATTERY REPLACEMENT

Battery replacement must only take place in a non-hazardous area using a 9V alkaline battery.

Simply slide off the battery cover on the back of the meter. Detach the battery from the terminals and attach a fresh 9V battery while paying attention to the correct polarity. Replace the battery and the cover.



ACCESSORIES

REAGENT SETS

HI 93748-01 Reagents for manganese LR (100 pcs of HI93748A-0, 20 mL of HI93748B-0, 100 mL of HI93748C-0, 20 mL of HI93703-51)

HI 93748-03 Reagents for manganese LR (300 pcs of HI93748A-0, 60 mL of HI93748B-0, 300 mL of HI93748C-0, 60 mL of HI93703-51)

OTHER ACCESSORIES

HI 710009 Blue rubber boot

HI 710010 Orange rubber boot

HI 721310 9V battery (10 pcs)

HI 731318 Tissue for wiping cuvets (4 pcs)

HI 731321 Glass cuvets (4 pcs)

HI 731325 Caps for cuvets (4 pcs)

HI 93703-50 Cuvets cleaning solution (230 mL).

CE DECLARATION OF CONFORMITY



Recommendations for Users

Before using these products, make sure that they are entirely suitable for the environment in which they are used.

Operation of these instruments in residential area could cause unacceptable interferences to radio and TV equipments, requiring the operator to take all necessary steps to correct interferences.

Any variation introduced by the user to the supplied equipment may degrade the instruments' EMC performance.

To avoid damages or burns, do not perform any measurement in microwave ovens.

PRINTED IN PORTUGAL

IST93748R1

06/98