



HI 190M-0

Mini Magnetic Stirrer with Speedsafe



PEWA
Messtechnik GmbH
Weidenweg 21
58239 Schwerte

Tel.: 02304-96109-0
Fax: 02304-96109-88
E-Mail: info@pewa.de
Homepage : www.pewa.de

- Compact size - The compact size of this stirrer allows users to maximize bench space for efficiency and safety.
- Safety - Speedsafe™ limits the maximum speed to 1000 rpm even if a load is suddenly removed.
- Built to last - The ABS housing of HI 190 M-0 resists most harmful chemicals in the lab.

The HI 190M-0 is compact and lightweight, so that lack of laboratory bench space is no longer a concern.

This stirrer incorporates electronic controls that allow the user to regulate the speed with greater precision. Often, in the lab, a sample is removed from the stirrer before reducing the speed. This would cause the motor of conventional equipment to accelerate until it is destroyed. This does not pose a problem with HANNA mini-stirrers. In fact, the internal Speedsafe™ mechanism will ensure that the maximum speed is never exceeded.

HI 190M-0 comes supplied with an ABS cover that will resist the harmful effects of chemicals that are accidentally spilled.

Order Information:

HI 190M-0 mini-stirrer with ABS cover, 12 Vdc is supplied with micro stir bar and instructions.

Specifications	Accessories
----------------	-------------

Maximum Stirring Capacity	1 liter (0.26 gallons)
Speed Range	Min. 100 rpm Max. 1000 rpm
Power Supply	12 Vdc
Installation Category	II
Cover Material	ABS plastic
Environment	0 to 50°C (32 to 122°F); RH max 95%
Dimensions	120 x 120 x 45 mm (4.8 x 4.8 x 1.8")
Weight	610 g (1.3 lb.)