

Weidenweg 21

58239 Schwerte

Tel.: 02304-96109-0 Fax: 02304-96109-88 E-Mail: info@pewa.de Homepage: www.pewa.de

HI 2209

pH/ORP Bench-Top Meter

The HI 2209 pH/mV Meter with Manual Temperature Compensation (MTC) provides a simple to use, cost effective method for measuring pH. Designed for durability and simplicity of use, this instrument features a large, easy to read LCD and user friendly dials.

In order to achieve maximum accuracy, the HI 2209 features a manual pH calibration at one or two points. Manual calibration enables the user to select the instrument's calibration points closer to the desired range of measurement, making them ideal for applications that require custom calibration points. (i.e.- In many applications, a standard calibration curve such as pH 7 or pH 4 is too far from the value of the sample to achieve the highest accuracy.)

The HI 2209 pH/mV Meter can also measure ion concentration (ISE) or ORP (oxidation reduction potential) in the extended mV range with optional electrodes.

FEATURES/BENEFITS:

- * Manual calibration for offset and slope
- * Manual Temperature Compensation (MTC)
- * Large, easy to read LCD
- * Built-in Solution Holders
- * Includes HI 1332B PEI body refillable probe

FEATURES IN DEPTH:

Manual pH calibration - This simple to use feature provides the ability to demonstrate the concept of offset and slope. It can be calibrated to any value within the measurement ranges and is less expensive than models with auto-calibration.

Manual Temperature Compensation (MTC) – MTC provides the ability to demonstrate the effect of temperature on pH measurement. It is simple to use and allows for different temperature corrections based on the sample being tested.

mV Range - Allows an ORP electrode to be used to measure REDOX reactions. The mV function can also be used with ISE electrodes for specific ion concentrations.

Large LCD – The new, larger LCD is bright and easy to read.

Built-in Solution Holders - The HI 2209 has solution holders built into the casing. This convenient feature saves space and prevents solutions from tipping over.

Recessed Connector Bay - Protects connectors from liquids.

PRIMARY APPLICATIONS:

- * Water Quality- Testing water quality
- * Education Teach students about pH

Electrodes

HI 1332B Use: general purpose. PEI body, double junction, refillable HI 3131B Glass body ORP electrode with platinum sensor, BNC connector and 1m (3.3') cable

Solutions

HI 5100-12 10.01 pH Buffer Solution (120 ml)
HI 5300-12 Electrode Storage Solution (120 ml)
HI 5400-12 4.01 pH Buffer Solution (120 ml)
HI 5700-12 7.01 pH Buffer Solution (120 ml)
HI7061M Electrode Cleaning Solution
HI 7082 Electrolyte Solution, 3.5M KCl, 4 x 30 mL bottle

Accessories

HI 76404 Electrode Holder HI 710005 Power Adapter, 115 Vac to 12 VDC, US Plug HI 710006 Power Adapter, 230 Vac to 12 VDC, European Plug HI 740036P 100 mL Plastic Beaker

Order Information:

 $\,$ HI 2209-01 (115V) and HI 2209-02 (230V) are supplied with HI 1332B pH electrode, 12 VDC adapter and instruction manual